

Mouse PRL (Prolactin) Antibody Pair (STJAP00540)

STJAP00540

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PRL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A & Antigen Affinity. Detection Purification: Protein A & Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	Prl
Uniprot ID	PRL_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse PRL protien. Detection Immunogen: Recombinant Mouse PRL protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	