

## Mouse RAN Kappa L (Tumor Necrosis Factor Ligand Superfamily Member 11) Antibody Pair (STJAP00541)

STJAP00541

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This RAN Kappa L antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Goat/Goat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">21943</a>
<b>Gene Symbol</b>	<a href="#">Tnfrsf11</a>
<b>Uniprot ID</b>	<a href="#">TNF11_MOUSE</a>
<b>Immunogen</b>	Capture immunogen: Recombinant Mouse RAN Kappa L protien. Detection Immunogen: Recombinant Mouse RAN Kappa L protien.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	