

Mouse TF (Tissue Factor) Antibody Pair (STJAP00546)

STJAP00546

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This TF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14066
Gene Symbol	F3
Uniprot ID	TF_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse TF protien. Detection Immunogen: Recombinant Mouse TF protien.
Immunogen Region	
Specificity	
Immunogen Sequence	