

Mouse bFGF/FGF2 (Fibroblast Growth Factor 2) Antibody Pair (STJAP00556)

STJAP00556

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This bFGF/FGF2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	
Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2247
Gene Symbol	FGF2
Uniprot ID	FGF2_HUMAN
Immunogen	Capture immunogen: Recombinant Mouse bFGF/FGF2 protien. Detection Immunogen: Recombinant Mouse bFGF/FGF2 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	