

## Mouse TRACP-5b (Tartrate-Resistant Acid Phosphatase Type 5) Antibody Pair (STJAP00558)

STJAP00558

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Antibody Pairs   |
| <b>Short Description</b> | This TRACP-5b antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| <b>Applications</b>      | ELISA  |
| <b>Host/Source</b>       | Rabbit/Rabbit  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                       |   |
|-----------------------|---|
| <b>Clonality</b>      |   |
| <b>Clone ID</b>       |   |
| <b>Concentration</b>  | 0.5 mg/mL   |
| <b>Conjugation</b>    | Unconjugated/Biotin   |
| <b>Purification</b>   | Capture Purification: Antigen Affinity.<br>Detection Purification: Antigen Affinity.  |
| <b>Dilution Range</b> |   |
| <b>Formulation</b>    | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5.<br>Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| <b>Isotype</b>        |   |
| <b>Storage</b>        |   |
| <b>Instruction</b>    | Store at -20°C for long term, and avoid freeze-thaw cycles.   |

### TARGET INFORMATION

|                    |  |
|--------------------|--|
| <b>Gene ID</b>     | 11433  |
| <b>Gene Symbol</b> | Acp5   |
| <b>Uniprot ID</b>  | PPA5_MOUSE   |
| <b>Immunogen</b>   | Capture immunogen: Recombinant Mouse TRACP-5b protien.<br>Detection Immunogen: Recombinant Mouse TRACP-5b protien. |
| <b>Immunogen</b>   |  |
| <b>Region</b>      |  |
| <b>Specificity</b> |  |
| <b>Immunogen</b>   |  |
| <b>Sequence</b>    |  |