

## Rat OC/BGP (Osteocalcin) Antibody Pair (STJAP00597)

STJAP00597

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This OC/BGP antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse/Mouse
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">25295</a>
<b>Gene Symbol</b>	<a href="#">Bglap</a>
<b>Uniprot ID</b>	<a href="#">OSTCN_RAT</a>
<b>Immunogen</b>	Capture immunogen: Recombinant Rat OC/BGP protien. Detection Immunogen: Recombinant Rat OC/BGP protien.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	