

Rat IgA (IgA) Antibody Pair (STJAP00599)

STJAP00599

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IgA antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Natural Rat IgA protien. Detection Immunogen: Natural Rat IgA protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	