

Rat DKK1 (DKK1) Antibody Pair (STJAP00607)

STJAP00607

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This DKK1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Sheep/Sheep
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Recombinant Rat DKK1 protien. Detection Immunogen: Recombinant Rat DKK1 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	