

Rat PAI1 (Plasminogen Activator Inhibitor 1) Antibody Pair (STJAP00614)

STJAP00614

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PAI1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24617
Gene Symbol	Serpine1
Uniprot ID	PAI1_RAT
Immunogen	Capture immunogen: Recombinant Rat PAI1 protien. Detection Immunogen: Recombinant Rat PAI1 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	