

Rat PTX-3/TSG-14 (PTX-3/TSG-14) Antibody Pair (STJAP00615)

STJAP00615

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PTX-3/TSG-14 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Recombinant Rat PTX-3/TSG-14 protien. Detection Immunogen: Recombinant Rat PTX-3/TSG-14 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	