

Rat IL-6R (Interleukin-6 Receptor Subunit Alpha) Antibody Pair (STJAP00622)

STJAP00622

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Antibody Pairs |
| Short Description | This IL-6R antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| Applications | ELISA |
| Host/Source | Goat/Goat |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity. |
| Dilution Range | |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| Isotype | |
| Storage | Store at -20°C for long term, and avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | |
| Gene Symbol | IL6r |
| Uniprot ID | IL6RA_RAT |
| Immunogen | Capture immunogen: Recombinant Rat IL-6R protein. Detection Immunogen: Recombinant Rat IL-6R protein. |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |