

## Human ACAT2 (ACAT2) CLIA Kit (STJC0001984)

STJC0001984

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This human ACAT2 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes in tissue homogenates, cell lysates and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Chemiluminescence
<b>Sensitivity</b>	18.75 pg/mL
<b>Detection Limit</b>	31.25-2000 pg/mL
<b>Assay Time</b>	3.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Instruction</b>	
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Tissue homogenates, cell lysates and other biological fluids
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	