

Human ADCY1 (ADCY1) CLIA Kit (STJC0001995)

STJC0001995

GENERAL INFORMATION

Product Type	Kits
Short Description	This human ADCY1 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes in tissue homogenates, cell lysates and other biological fluids.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Chemiluminescence
Sensitivity	0.09 ng/mL
Detection Limit	0.16-10 ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Tissue homogenates, cell lysates and other biological fluids
Gene	
Symbol	
Uniprot ID	