

**Mouse IL-2sR Beta ;CD122 (IL-2sR Beta ;CD122) CLIA Kit (STJC0002034)**

STJC0002034

**GENERAL INFORMATION**

<b>Product Type</b>	Kits
<b>Description</b>	This mouse IL-2sR Beta ;CD122 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes in serum, plasma and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Detection</b>	Chemiluminescence
<b>Sensitivity</b>	9.38 pg/mL
<b>Detection Limit</b>	15.63-1000 pg/mL
<b>Assay Time</b>	3.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Instruction</b>	
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

**TARGET INFORMATION**

<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Gene</b>	
<b>Symbol</b>	
<b>Uniprot ID</b>	