

Human IFNL2 (IFNL2) Sandwich ELISA Kit Kit (STJE0000082)

STJE0000082

GENERAL INFORMATION

Product Type	Kits
Short Description	
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	7.8 pg/mL
Detection Limit	156-1000 pg/mL
Assay Time	4 hours
Storage	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and cell culture supernates
Gene Symbol	IFNL2
Uniprot ID	IFNL2_HUMAN