

Human FNDC5 (FNDC5) Sandwich ELISA Kit (STJE0000767)

STJE0000767

GENERAL INFORMATION


| | |
|--------------------------|-------|
| Product Type | Kits |
| Short Description | |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 6.9 pg/mL |
| Detection Limit | 312-2000 pg/ml |
| Assay Time | 4 hours |
| Storage | Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month. |
| Instruction Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma and cell culture supernates |
| Gene Symbol | FNDC5 |
| Uniprot ID | FNDC5_HUMAN |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081