

Human HIF1a (HIF1a) Sandwich ELISA Kit (STJE0000920)

STJE0000920

GENERAL INFORMATION


Product Type	Kits
Short Description	
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	35.1 pg/mL
Detection Limit	78-5000 pg/ml
Assay Time	4 hours
Storage	Unopened kits can be stored at 2-8°C for 1 year, and opened products must be used within 1 month.
Instruction Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and cell culture supernates
Gene Symbol	HIF1A
Uniprot ID	HIF1A_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081