

Human PROS (Protein S) Sandwich ELISA Kit (STJE0002765)

STJE0002765

GENERAL INFORMATION


Product Type	Kits
Short Description	This PROS Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human plasma.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.19ng/mL
Detection Limit	0.31~20ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Plasma
Gene Symbol	PROS1
Uniprot ID	PROS_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081