

## Human IL27 (Interleukin-27 subunit alpha) Sandwich ELISA Kit (STJE0003630)

STJE0003630

### GENERAL INFORMATION


|                          |                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                         |
| <b>Short Description</b> | This IL27 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids. |
| <b>Applications</b>      | ELISA                                                                                                                                                        |
| <b>Reactivity</b>        | Human                                                                                                                                                        |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | Colormetric                                                                                                                                                                                |
| <b>Sensitivity</b>     | 18.75pg/mL                                                                                                                                                                                 |
| <b>Detection Limit</b> | 31.25~2000pg/mL                                                                                                                                                                            |
| <b>Assay Time</b>      | 3.5h                                                                                                                                                                                       |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |                                                                                                                                                                                            |
| <b>Note</b>            | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.                                                                          |

### TARGET INFORMATION

|                     |                                           |
|---------------------|-------------------------------------------|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  | <a href="#">IL27</a>                      |
| <b>Uniprot ID</b>   | <a href="#">IL27A_HUMAN</a>               |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081