

## Mouse RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) Sandwich ELISA Kit (STJE0004209)

STJE0004209

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This RANTES Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	18.75pg/mL
<b>Detection Limit</b>	31.25~2000pg/mL
<b>Assay Time</b>	3.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Gene Symbol</b>	<a href="#">Ccl5</a>
<b>Uniprot ID</b>	<a href="#">CCL5_MOUSE</a>