

Mouse GRO Gamma/CXCL3 (Growth Regulated Oncogene Gamma) Sandwich ELISA Kit (STJE0004293)

STJE0004293

GENERAL INFORMATION

Product Type	Kits
Short Description	This GRO Gamma/CXCL3 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	4,69pg/mL
Detection Limit	7.81~500pg/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Cxcl3
Uniprot ID	CXCL3_MOUSE