

Mouse Cxcl13 (C-X-C motif chemokine 13) Sandwich ELISA Kit (STJE0004328)

STJE0004328

GENERAL INFORMATION

Product Type	Kits
Short Description	This Cxcl13 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	9.38pg/mL
Detection Limit	15.63~1000pg/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Cxcl13
Uniprot ID	CXL13_MOUSE