

**Mouse Ctla4 (Cytotoxic T-lymphocyte protein 4) Sandwich ELISA Kit (STJE0004495)**

STJE0004495

**GENERAL INFORMATION**

<b>Product Type</b>	Kits
<b>Description</b>	This Ctla4 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	75pg/mL
<b>Detection Limit</b>	125~8000pg/mL
<b>Assay Time</b>	3.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

**TARGET INFORMATION**

<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Gene Symbol</b>	<a href="#">Ctla4</a>
<b>Uniprot ID</b>	<a href="#">CTLA4_MOUSE</a>

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081