

Mouse Ccl7 (C-C motif chemokine 7) Sandwich ELISA Kit (STJE0004785)

STJE0004785

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This Ccl7 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids. |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | Colormetric |
| Sensitivity | 9.38pg/mL |
| Detection Limit | 15.63~1000pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene Symbol | Ccl7 |
| Uniprot ID | CCL7_MOUSE |