

**Mouse MMP-9 (Matrix Metalloproteinase 9) Sandwich ELISA Kit (STJE0005368)**

STJE0005368

**GENERAL INFORMATION**

|                     |   |
|---------------------|---|
| <b>Product Type</b> | Kits  |
| <b>Description</b>  | This MMP-9 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma and other biological fluids. |
| <b>Applications</b> | ELISA   |
| <b>Reactivity</b>   | Mouse   |

**PRODUCT PROPERTIES**

|                        |  |
|------------------------|--|
| <b>Detection</b>       | Colormetric  |
| <b>Sensitivity</b>     | 0.94ng/mL  |
| <b>Detection Limit</b> | 1.56~100ng/mL  |
| <b>Assay Time</b>      | 3.5h   |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Note</b>            | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.  |

**TARGET INFORMATION**

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  | <a href="#">Mmp9</a>                      |
| <b>Uniprot ID</b>   | <a href="#">MMP9_MOUSE</a>                |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081