

Rat Uchl1 (Ubiquitin carboxyl-terminal hydrolase isozyme L1) Sandwich ELISA Kit (STJE0006430)

STJE0006430

GENERAL INFORMATION


Product Type	Kits
Short Description	This Uchl1 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	37.5pg/mL
Detection Limit	62.5~4000pg/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Uchl1
Uniprot ID	UCHL1_RAT



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081