

Rat Havcr1 (Hepatitis A virus cellular receptor 1 homolog) Sandwich ELISA Kit (STJE0006554)

STJE0006554

GENERAL INFORMATION

Product Type	Kits
Short Description	This Havcr1 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.19ng/mL
Detection Limit	0.31~20ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	Havcr1
Uniprot ID	HAVR1_RAT