

8-epi-PGF2 Alpha (8-Epi-Prostaglandin F2 Alpha) Competitive ELISA Kit (STJE0006589)

STJE0006589

GENERAL INFORMATION


Product Type	Kits
Short Description	This 8-epi-PGF2 Alpha Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in all serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	All

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	9.38pg/mL
Detection Limit	15.63~1000pg/mL
Assay Time	2.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081