

## Human CD8 (Cluster of Differentiation 8) Sandwich ELISA Kit (STJE0006653)

STJE0006653

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This CD8 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	0.75ng/mL
<b>Detection Limit</b>	1.25~80ng/mL
<b>Assay Time</b>	3.5h
<b>Storage</b>	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Instruction</b>	
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Gene Symbol</b>	<a href="#">CD8B</a> <a href="#">CD8A</a>
<b>Uniprot ID</b>	<a href="#">CD8B_HUMAN</a> <a href="#">CD8A_HUMAN</a>