

Sars-Cov-2 SARS-CoV-2 (SARS-CoV-2 Nucleocapsid Protein) Sandwich ELISA Kit (STJE0006698)

STJE0006698

GENERAL INFORMATION

Product Type	Kits
Short Description	This SARS-CoV-2 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.12ng/mL
Detection Limit	0.39~25ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	