

Sars-Cov-2 SARS-CoV-2 Spike Protein Total Antibody (SARS-CoV-2 Spike Protein Total Antibody) Sandwich ELISA Kit (STJE0006703)

STJE0006703

GENERAL INFORMATION

Product Type	Kits
Short Description	This SARS-CoV-2 Spike Protein Total Antibody Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in sars-cov-2 serum and plasma.
Applications	ELISA
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	
Detection Limit	
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	
Uniprot ID	