

## Human LCN2 (Neutrophil gelatinase-associated lipocalin) Sandwich ELISA Kit (STJE0006705)

STJE0006705

### GENERAL INFORMATION

|                          |                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                         |
| <b>Short Description</b> | This LCN2 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids. |
| <b>Applications</b>      | ELISA                                                                                                                                                        |
| <b>Reactivity</b>        | Human                                                                                                                                                        |

### PRODUCT PROPERTIES

|                            |                                                                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>           | Colormetric                                                                                                                                                                                |
| <b>Sensitivity</b>         | 18.75pg/ml                                                                                                                                                                                 |
| <b>Detection Limit</b>     | 31.25~2000 pg/mL                                                                                                                                                                           |
| <b>Assay Time</b>          | 3.5h                                                                                                                                                                                       |
| <b>Storage Instruction</b> | If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Note</b>                | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.                                                                          |

### TARGET INFORMATION

|                     |                                           |
|---------------------|-------------------------------------------|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  | <a href="#">LCN2</a>                      |
| <b>Uniprot ID</b>   | <a href="#">NGAL_HUMAN</a>                |