

## Bovine Cortisol (Cortisol) Competitive ELISA Kit (STJE0006716)

STJE0006716

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This Cortisol Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in bovine serum and plasma. |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Bovine  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Detection</b>           | Colorimetric   |
| <b>Sensitivity</b>         | 2.54ng/mL  |
| <b>Detection Limit</b>     | 6.25~400ng/mL  |
| <b>Assay Time</b>          | 1.5h   |
| <b>Storage Instruction</b> | If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual. |
| <b>Note</b>                | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.  |

### TARGET INFORMATION

|                     |                  |
|---------------------|------------------|
| <b>Sample Types</b> | Serum and plasma |
| <b>Gene Symbol</b>  |                  |
| <b>Uniprot ID</b>   |                  |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081