

Horse Horse E2 (Estradiol) Competitive ELISA Kit (STJE0006723)

STJE0006723

GENERAL INFORMATION

Product Type	Kits
Short Description	This Horse E2 Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in horse serum and plasma.
Applications	ELISA
Reactivity	Horse

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	3.46pg/mL
Detection Limit	7.81~500pg/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	
Uniprot ID	