

## Mouse CORT (Corticosterone) Competitive ELISA Kit (STJE0006724)

STJE0006724

### GENERAL INFORMATION


<b>Product Type</b>	Kits
<b>Short Description</b>	This CORT Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma, urine and saliva.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	Colorimetric
<b>Sensitivity</b>	4.98ng/mL
<b>Detection Limit</b>	12.5~800ng/mL
<b>Assay Time</b>	1.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Serum, Plasma, Urine and Saliva
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081