

Rabbit Cortisol (Cortisol) Competitive ELISA Kit (STJE0006733)

STJE0006733

GENERAL INFORMATION


Product Type	Kits
Short Description	This Cortisol Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rabbit serum and plasma.
Applications	ELISA
Reactivity	Rabbit

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	2.87ng/mL
Detection Limit	6.25~400ng/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081