

Human IL6 (Interleukin-6) Sandwich ELISA Kit (STJE0006744)

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL6 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum and plasma.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	1.02pg/mL
Detection Limit	3.13-200pg/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN