

Human CRP (C-reactive protein) Sandwich ELISA Kit (STJE0006745)

STJE0006745

GENERAL INFORMATION

Product Type	Kits
Short Description	This CRP Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum and plasma.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	52.1pg/mL
Detection Limit	125-8000pg/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	CRP
Uniprot ID	CRP_HUMAN