

Human FE (Ferritin) Sandwich ELISA Kit (STJE0006746)

STJE0006746

GENERAL INFORMATION

Product Type	Kits
Short Description	This FE Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum and plasma.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	24.31pg/mL
Detection Limit	62.5-4000pg/mL
Assay Time	1.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	FTL
	FTH1
Uniprot ID	FRIL_HUMAN
	FRIH_HUMAN