

## Human CXCL8 (Interleukin-8) Sandwich ELISA Kit (STJE0006755)

STJE0006755

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This CXCL8 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum and plasma.
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Colormetric
<b>Sensitivity</b>	9.38pg/mL
<b>Detection Limit</b>	15.63-1000pg/mL
<b>Assay Time</b>	1.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
<b>Note</b>	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

### TARGET INFORMATION

<b>Sample Types</b>	Serum and plasma
<b>Gene Symbol</b>	<a href="#">CXCL8</a>
<b>Uniprot ID</b>	<a href="#">IL8_HUMAN</a>