

Human CEACAM5 (Carcinoembryonic antigen-related cell adhesion molecule 5) Sandwich ELISA Kit (STJE0006757)

STJE0006757

GENERAL INFORMATION

Product Type	Kits
Short Description	This CEACAM5 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum and plasma.
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.47ng/mL
Detection Limit	0.78-50ng/mL
Assay Time	1.5h
Storage	If unopened the kit may be stored at 2-8°C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Instruction	
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

TARGET INFORMATION

Sample Types	Serum and plasma
Gene Symbol	CEACAM5
Uniprot ID	CEAM5_HUMAN