

## Mouse ADCY10 (Adenylate Cyclase 10, Soluble) Sandwich ELISA Kit (STJE0006911)

STJE0006911

### GENERAL INFORMATION


|                          |                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                    |
| <b>Short Description</b> | This ADCY10 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse tissue homogenates or other biological fluids.. |
| <b>Applications</b>      | ELISA                                                                                                                                                                   |
| <b>Reactivity</b>        | Mouse                                                                                                                                                                   |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | 0.312-20ng/mL                                                                                                                                                    |
| <b>Sensitivity</b>     | 0.132ng/mL                                                                                                                                                       |
| <b>Detection Limit</b> | 0.312-20ng/mL                                                                                                                                                    |
| <b>Assay Time</b>      | 4.5 hrs                                                                                                                                                          |
| <b>Storage</b>         | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Instruction</b>     |                                                                                                                                                                  |
| <b>Note</b>            | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                          |

### TARGET INFORMATION

|                     |                                                |
|---------------------|------------------------------------------------|
| <b>Sample Types</b> | tissue homogenates or other biological fluids. |
| <b>Gene Symbol</b>  | <a href="#">Adcy10</a>                         |
| <b>Uniprot ID</b>   | <a href="#">ADCYA_MOUSE</a>                    |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081