

Mouse PF4 (Platelet Factor 4) Sandwich ELISA Kit Ready-To-Use Kit (STJE0007025)

STJE0007025

GENERAL INFORMATION


Product Type	Kits
Short Description	This PF4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, platelet-poor plasma, cell culture supernates or other biological fluids..
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	27pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, platelet-poor plasma, cell culture supernates or other biological fluids.
Gene	Pf4
Symbol	
Uniprot ID	PLF4_MOUSE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081