

Rat FGg (Fibrinogen Gamma) Sandwich ELISA Kit (STJE0008923)

STJE0008923

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This FGg Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat plasma.. |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | 78.1-5000ng/mL |
| Sensitivity | 27ng/mL |
| Detection Limit | 78.1-5000ng/mL |
| Assay Time | 4.5 hrs |
| Storage | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Instruction | table provided with the kit. |
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |

TARGET INFORMATION

| | |
|---------------------|--------------------------|
| Sample Types | plasma. |
| Gene Symbol | Fgg |
| Uniprot ID | FIBG_RAT |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081