

Human HBg1 (Hemoglobin Gamma 1) Competitive ELISA Kit (STJE0009399)

STJE0009399

GENERAL INFORMATION


Product Type	Kits
Short Description	This HBg1 Competitive ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma or erythrocyte lysates..
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	468-30000ng/mL
Sensitivity	122.6ng/mL
Detection Limit	468-30000ng/mL
Assay Time	2.5 hrs
Storage Instruction	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

TARGET INFORMATION

Sample Types	serum, plasma or erythrocyte lysates.
Gene Symbol	HBG1
Uniprot ID	HBG1_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081