

## Mouse ITGa4 (Integrin Alpha 4) Sandwich ELISA Kit (STJE0010082)

STJE0010082

### GENERAL INFORMATION


|                          |                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                                 |
| <b>Short Description</b> | This ITGa4 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse tissue homogenates, cell lysates or other biological fluids.. |
| <b>Applications</b>      | ELISA                                                                                                                                                                                |
| <b>Reactivity</b>        | Mouse                                                                                                                                                                                |

### PRODUCT PROPERTIES

|                            |                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>           | 0.156-10ng/mL                                                                                                                                                    |
| <b>Sensitivity</b>         | 0.071ng/mL                                                                                                                                                       |
| <b>Detection Limit</b>     | 0.156-10ng/mL                                                                                                                                                    |
| <b>Assay Time</b>          | 4.5 hrs                                                                                                                                                          |
| <b>Storage Instruction</b> | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                          |

### TARGET INFORMATION

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| <b>Sample Types</b> | tissue homogenates, cell lysates or other biological fluids. |
| <b>Gene Symbol</b>  | <a href="#">Itga4</a>                                        |
| <b>Uniprot ID</b>   | <a href="#">ITA4_MOUSE</a>                                   |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081