

## Rat PYGL (Glycogen Phosphorylase, Liver) Sandwich ELISA Kit (STJE0011788)

STJE0011788

### GENERAL INFORMATION


|                          |                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                               |
| <b>Short Description</b> | This PYGL Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat serum, plasma, tissue homogenates or other biological fluids.. |
| <b>Applications</b>      | ELISA                                                                                                                                                                              |
| <b>Reactivity</b>        | Rat                                                                                                                                                                                |

### PRODUCT PROPERTIES

|                            |                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>           | 1.56-100ng/mL                                                                                                                                                    |
| <b>Sensitivity</b>         | 0.56ng/mL                                                                                                                                                        |
| <b>Detection Limit</b>     | 1.56-100ng/mL                                                                                                                                                    |
| <b>Assay Time</b>          | 4.5 hrs                                                                                                                                                          |
| <b>Storage Instruction</b> | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                          |

### TARGET INFORMATION

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| <b>Sample Types</b> | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Gene Symbol</b>  | <a href="#">Pygl</a>                                          |
| <b>Uniprot ID</b>   | <a href="#">PYGL_RAT</a>                                      |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081