

Mouse Dyn (Big Dynorphin) Competitive ELISA Kit Ready-To-Use Kit (STJE0014915)

STJE0014915

GENERAL INFORMATION


| | |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product | Kits |
| Type | |
| Short Description | This Dyn Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma, tissue homogenates or other biological fluids.. |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | 4.68-300pg/mL |
| Sensitivity | 2.4pg/mL |
| Detection Limit | 4.68-300pg/mL |
| Assay Time | 2.5 hrs |
| Storage | The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Instruction | |
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |

TARGET INFORMATION

| | |
|-------------------|---------------------------------------------------------------|
| Sample | serum, plasma, tissue homogenates or other biological fluids. |
| Types | |
| Gene | Pdyn |
| Symbol | |
| Uniprot ID | PDYN_MOUSE |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081