

## Guinea Pig GT (pig Gastrin) Sandwich ELISA Kit Ready-To-Use Kit (STJE0015734)

STJE0015734

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product</b>           | Kits  |
| <b>Type</b>              |   |
| <b>Short Description</b> | This GT Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in guinea pig serum, plasma or other biological fluids.. |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Guinea Pig  |

### PRODUCT PROPERTIES

|                        |  |
|------------------------|--|
| <b>Detection</b>       | 15.6-1000pg/mL   |
| <b>Sensitivity</b>     | 5.9pg/mL   |
| <b>Detection Limit</b> | 15.6-1000pg/mL   |
| <b>Assay Time</b>      | 3 hrs  |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Instruction</b>     |  |
| <b>Note</b>            | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.  |

### TARGET INFORMATION

|                   |   |
|-------------------|---|
| <b>Sample</b>     | serum, plasma or other biological fluids. |
| <b>Types</b>      |   |
| <b>Gene</b>       |   |
| <b>Symbol</b>     |   |
| <b>Uniprot ID</b> |   |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081