

**Bovine IL12A (Interleukin 12A) Sandwich ELISA Kit Ready-To-Use Kit (STJE0016143)**

STJE0016143

**GENERAL INFORMATION**

<b>Product</b>	Kits
<b>Type</b>	
<b>Short</b>	This IL12A Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in
<b>Description</b>	bovine serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids..
<b>Applications</b>	ELISA
<b>Reactivity</b>	Bovine

**PRODUCT PROPERTIES**

<b>Detection</b>	15.625-1000pg/mL
<b>Sensitivity</b>	5.9pg/mL
<b>Detection</b>	15.625-1000pg/mL
<b>Limit</b>	
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 12 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6
<b>Instruction</b>	months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

**TARGET INFORMATION**

<b>Sample</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Types</b>	
<b>Gene</b>	
<b>Symbol</b>	
<b>Uniprot ID</b>	